


**Fotona**  
Ultra Performance Lasers™

**LIGHTWALKER™**  
Dual Wavelength Nd:YAG & Er:YAG Laser




**SWEEPS®**  
Endo Laser Treatment

SWEEPS® is a revolutionary method for cleaning the root canal system, utilizing the power of the Er:YAG laser to create powerful, deeply penetrating shockwaves within the cleaning and debriding solutions in the root canal.



**BEFORE**  
Bacteria left in  
Traditional Endo



**AFTER**  
LightWalker Laser  
Endo Treatment

- Improved debridement
- Precise deep penetration
- Superior activation of irrigants
- Minimally invasive
- Easily reach accessory & lateral canals
- Locate difficult to find canals



Before my LightWalker, there is no way I truly cleaned canal systems. I cleaned canals but there was no way I was able to get into dentinal tubules, areas where bacteria was hiding.



Scan to watch Dr. Haas' webinar on SWEEPS Endo

**- Manor Haas, DDS**  
Toronto, Canada

"Traditional root canal treatments do not clean everything out. They only clean a small, tiny percentage of the dentinal tubules, leaving a ton of bacteria and toxins behind... but with the laser, you're sterilizing cell layer by cell layer, and it's really incredible."



**- Valerie Kanter, DDS**  
Integrative Endodontics, Beverly Hills, CA

Scan the QR code to learn more about how Endodontists are using SWEEPS and how the LightWalker can benefit your patients and your practice!